



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybierański  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Medical Stool Coburg Medicalift 22011 (New)

- New,
- Made in Germany,
- Multi-purpose stool suitable for surgery and outpatient departments,
- Safety base made of stainless steel with 5 anti-static twin castors for hard floors,
- Chrome-plated foot ring,
- Height adjustment using a safety gas spring controlled by means of stainless steel foot-controlled disc,
- Seat height: 60 - 83 cm,
- Standard version: black saddle seat,
- **Additional Equipment (extra paid):**
  - backrest mod. 2\*,
  - rounded seat Ø35cm\*,
  - rounded seat Ø38cm\*,
  - moulded seat 1103\*,
  - comfortable upholstery as per our sample card,
  - special colours subject to a surcharge,
- Technical Passport,
- 24-month warranty

