



Medical Stool Coburg Medicalift 22012 (New)

- New,
- Made in Germany,
- Multi-purpose stool suitable for surgery, anaesthesia and outpatient departments,
- Stainless steel security base with 5 twin castors for hard floors,
- Height adjustment using a safety gas spring controlled by means of a stainless steel foot-controlled disc,
- Rounded seat with anti-static black upholstery,
- Seat Height:
 - H200: 47 - 66 cm,
 - H250: 52 - 76 cm,
- **Additional Equipment (extra paid):**
 - back rest with anti-static black upholstery,
 - special colours subject to a surcharge,
 - central circular stop bar operated by placing your feet on the footplate,
- Technical Passport,
- 24-month warranty

