



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Coburg Exaflex 6115 operating table (New)**

- A brand new product,
- Made in Germany,
- **The basic version of the table has:**
  - A headrest (standard),
  - A dividable footrest (standard),
- Portable multifunctional surgical table for general use in hospitals and ambulatory care units,
- Also suitable for military use due to its relatively light weight and ease of disassembly. The table can be disassembled and packed into two aluminum cases,
- The table is equipped as standard with chrome-nickel steel rails located on both sides of the seat, backrest and footrest,
- The table is covered with anti-static, washable upholstery which is resistant to all kinds of disinfectants, and does not absorb any liquids, including blood,
- Height regulation by a hydraulic pump lever,
- Backrest regulation by the use of gas spring,
- Lateral tilting of the table and Trendelenburg and anti-Trendelenburg positions by the use of a crank mechanism,
- Table plates made of radio transparent polyamide material,
- Movable base made of stainless steel with four double anti-static wheels (Ø100mm ) and central brake,
- The table has the possibility of mounting an X-ray cassette,
- **Additional options (for extra fee):**
  - Splittable footrest, foldable in the knee area,
  - Patient hand support base,
  - Equalization connection wire,
- **Additional accessories for the table, visible in the photo (for extra fee),**
  - **24 months for the international market,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
  - **24 months for the domestic market (Poland),**
- Financing possibilities (Poland only): Installments, leasing, loan,

