



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NGP 400 P vertical negatoscope (New)

- A brand new product,
- Made in Poland,
- The vertical negatoscope is one type of negatoscopes used in top-quality medical equipment,
- The negatoscope is successfully used for accurate and precise review of X-ray images,
- **Advantages of using the vertical NGP 400 P negatoscope:**
 - Emits daylight (Day Light),
 - Evenly distributes light across the screen surface,
 - Easy and fast service, including replacement of light sources,
 - Attractive and durable case,
 - It has all required certifications, including entry in the Register of Medicinal Products, Medical Devices and Biocidal Products, CE declaration of conformity, ISO 13485:2003, ISO 9001:2008,
 - The roll holder strip can be made in any color,
- **Technical data:**
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 345 VA,
 - Screen luminance: 4100 cd/m²±25%,
 - Luminance regulation: optional,
 - Screen dimensions: 36 x 17 cm,
 - Number of frames: 1,
 - Protection class: I,
 - Dimensions: 410 x 1270 x 112 mm,
 - Weight: 19 kg,
- Provided with an up-to-date Technical Passport,
- Guarantee: 12 months,

