



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NGP 600 P vertical negatoscope (New)

- A brand new product,
- Made in Poland,
- The negatoscope provides an excellent X-ray image that, thanks to its precision, shows even the smallest lesions in soft tissues and any micro degeneration of other kinds,
- **Advantages of using the vertical NGP 600 P negatoscope:**
 - They have all the necessary certificates,
 - The roll holder strip can be made in any color,
 - Emits daylight with a color temperature of 6500°K,
 - Evenly distributes light over the entire screen surface,
 - It has easy and fast service, including light source replacement,
 - It is made in an attractive and durable case,
- **Technical data:**
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 550 VA,
 - Screen luminance: 4100 cd/m²±25%,
 - Screen dimensions: 72 x 87 cm,
 - Number of frames: 1,
 - Protection class: I,
 - Dimensions: 760 x 965 x 112 mm,
 - Weight: 32 kg,
- Provided with an up-to-date Technical Passport,
- Guarantee: 12 months,

