



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NGP 30 hf negatoscope (New)

- A brand new product,
- Made in Poland,
- The NGP 30 HF model is a negatoscope that allows for very accurate and convenient review of X-ray images,
- Allows to see even the smallest soft tissue lesions and any micro degeneration of other types.
- **Advantages of NGP 30 hf negatoscope:**
 - This model emits daylight with a color temperature of 6500K,
 - NGP series devices feature excellent uniformity of screen illumination,
 - Medical certifications,
 - Type of light emitted: daylight with a color temperature of 6500K,
 - The case is powder coated, which makes it easy to keep the device clean
- **Technical data:**
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 195 W,
 - Screen luminance: 6000 cd/m²±25%,
 - Luminance regulation: No,
 - Screen dimensions: 108 x 43 cm,
 - Number of frames: 3,
 - Protection class: I,
 - Dimensions: 1100 x 510 x 92 mm,
 - Weight 14,1 kg,
- Provided with an up-to-date Technical Passport,
- Guarantee: 12 months,

