



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NGP 41 hf negatoscope (New)

- A brand new product,
- Made in Poland,
- The NGP 41 HF model has smooth luminescence, which allows the person using the device not to feel the sensation of light shaking,
- The negatoscope has all the necessary certifications, including registration in the Register of Medicinal Products, Medical Devices and Biocidal Products, ISO 13485:2003, ISO 9001:2008 and a CE Declaration of Conformity, which confirm its superior quality,
- **Advantages of NGP 41 hf negatoscope:**
 - Screen division: 4 frames for viewing X-ray negatives,
 - Fashionable design of the device,
 - Easy installation,
 - Highest quality proven by medical certificates,
 - Screen evenly illuminated at every point,
 - Powder coating of the case allows the device to be kept clean,
 - Emits daylight with a color temperature of 6500°K,
- **Technical data:**
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 260 W,
 - Screen luminance: 6000 cd/m²±25%,
 - Luminance regulation: Yes,
 - Screen dimensions: 144 x 43 cm,
 - Number of frames: 4,
 - Protection class: I,
 - Dimensions: 1460 x 510 x 92 mm,
 - Weight: 19,6 kg,
- Provided with an up-to-date Technical Passport,
- Guarantee: 12 months,

