



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NGP 400 Z negatoscope (New)

- A brand new product,
- Made in Poland,
- **Advantages of using the NGP 400 Z negatoscope:**
 - Four blinds allow easy "blending" of the negatoscope screen making it easier to describe images without side glare,
 - The device also allows "blending" near the edges on small areas of images,
 - Device in high-frequency version,
 - Very fast ignition of fluorescent lamps,
 - Negatoscope does not generate stroboscopic flicker effect, thus protecting radiologist's eyes,
 - Saves electricity by more than 25%,
 - The ability to regulate the intensity and the use of blinds facilitate accurate description of details such as micro calcifications, this is important, for example, when describing mammography,
 - Uniform illumination of the screen,
 - Negatoscope emits day light (Day Light) 6500K,
 - A handy picture holder is mounted around the perimeter of the screen,
 - The slim, compact case of the device is protected with a powder coating that is resistant to commonly used disinfectants,
 - Negatoscopes have all required certificates,
- **Technical datane:**
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 260 VA,
 - Screen luminance: >4500 cd / m²,
 - Luminance control: YES,
 - Screen dimensions: 140 x 43cm,
 - Number of frames: 1,
 - Protection class: I,
 - Dimensions: 1740 x 690 x 120 mm,
 - Weight: 57 kg,
- Provided with an up-to-date Technical Passport,
- Guarantee: 12 months,

