



NBV 30 N germicidal lamp (New)

- A brand new product,
- Made in Poland,
- NBV 30 N direct-acting germicidal lamp without a timer,
- This device emits UV-C radiation with a wavelength of 253.7 nm,
- The lamp irreversibly inactivates bacteria, viruses, molds, fungi and all other microorganisms,
- **The use of the device:**
 - Operating rooms,
 - Treatment and bandaging rooms,
 - Patient rooms,
 - Clinics,
 - Laboratories,
 - Pharmacies,
 - Pharmaceutical industry,
 - Food industry,
 - Cosmetics industry,
- **Advantages of NBV 30 N germicidal lamp:**
 - Type: N - universal, wall and ceiling direct acting germicidal lamp,
 - Purpose: to raise and maintain the level of microbiological cleanliness in rooms of various types,
 - Purpose: irreversible deactivation of all kinds of bacteria, viruses, molds, fungi and microorganisms,
 - Effect: effective disinfection of air and exposed surface,
 - Purpose: the use in operating rooms, treatment and bandaging rooms, operating rooms, clinics, laboratories and pharmacies, as well as in the pharmaceutical, food and cosmetic industries,
 - Additional information: this lamp is equipped with a high-quality aluminum reflector; in addition, it has all the required certificates, and its effectiveness is confirmed by tests,
- **Technical data:**
 - Wall mounted lamp,
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 33 W,
 - UV-C emitting element: 1 x 30 W,
 - UV-C irradiance at 1 m distance: 2.3 W / m²,
 - Radiant life: 9000 h,
 - Disinfectable area: 12 - 15 m²,
 - Lamp rotation (exposure angle can be adjusted): 200°,
 - Anti-schock safety class: I,
 - Case type: IP 20,
 - Operation type: continuous,
 - Chamber's dimensions: 960 x 85 x 135 mm,
 - Chamber's weight: 3 kg,
 - The extension's length: 120 mm,
- It has a valid Technical Passport,
- Guarantee: 12 months,

