



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **NBV 30 S germicidal lamp (New)**

- A brand new product,
- Made in Poland,
- NBV 30 S direct-acting germicidal lamp without a timer,
- This device emits UV-C radiation with a wavelength of 253.7 nm,
- The lamp irreversibly inactivates bacteria, viruses, molds, fungi and all other microorganisms,
- **The use of the device:**
  - Operating rooms,
  - Treatment and bandaging rooms,
  - Patient rooms,
  - Clinics,
  - Laboratories,
  - Pharmacies,
  - Pharmaceutical industry,
  - Food industry,
  - Cosmetics industry,
- **Advantages of NBV 30 S germicidal lamp:**
  - Mounting type: C - ceiling,
  - With this germicidal lamp, you will raise the level of microbiological purity,
  - Emitted radiation: UV-C with a wavelength of 253.7 nm,
  - When using this lamp, viruses, molds, fungi and microorganisms are irreversibly deactivated,
  - Purpose: effective disinfection of air and surfaces,
  - Reflector: this element is made of aulumin, which guarantees its highest quality,
  - Purpose: this germicidal lamp is used in medical institutes, hospitals and clinics, as well as in pharmaceutical, food and cosmetic industries,
  - The effectiveness of this lamp has been confirmed by tests; it is also confirmed by numerous certificates,
- **Technical data:**
  - Ceiling mounted lamp,
  - Supply voltage: 230 V 50 Hz,
  - Power consumption: 33 VA,
  - UV-C emitting element: 1 x 30 W,
  - UV-C irradiance at 1 m distance: 2.3 W / m<sup>2</sup>,
  - Radiant life: 9000 h,
  - Disinfectable area: 12 - 15 m<sup>2</sup>,
  - Lamp rotation (exposure angle can be adjusted): 230°,
  - Anti-schock safety class: I,
  - Case type: IP 20,
  - Operation type: continuous,
  - Chamber's dimensions: 960 x 85 x 135 mm,
  - Chamber's weight: 3 kg,
  - The extension's length: 500 mm,
- It has a valid Technical Passport,
- Guarantee: 12 months,

