



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NBV 30 P germicidal lamp (New)

- A brand new product,
- Made in Poland,
- NBV 30 P direct-acting germicidal lamp without a timer,
- This device emits UV-C radiation with a wavelength of 253.7 nm,
- The lamp irreversibly inactivates bacteria, viruses, molds, fungi and all other microorganisms,
- **The use of the device:**
 - Operating rooms,
 - Treatment and bandaging rooms,
 - Patient rooms,
 - Clinics,
 - Laboratories,
 - Pharmacies,
 - Pharmaceutical industry,
 - Food industry,
 - Cosmetics industry,
- **Advantages of NBV 30 P germicidal lamp:**
 - Mounting type: M – mobile, on a tripod,
 - The device has passed tests that confirm its superior quality,
 - These lamps emit UV-C radiation with a wavelength of 253.7 nm,
 - Bacteria, molds, fungi and viruses exposed to this lamp are deactivated,
 - The most important task of this germicidal lamp is to disinfect air and surfaces,
 - Albumen reflector proves the high quality of the device,
 - The lamp also has all the required certificates,
 - This device is designed for the use in operating rooms, treatment and bandaging rooms, patient rooms, clinics, laboratories, pharmaceutical, food and cosmetic industries,
- **Technical data:**
 - The lamp is placed on a mobile tripod,
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 33 W,
 - UV-C emitting element: 1 x 30 W,
 - UV-C irradiance at 1 m distance: 2.3 W / m²,
 - Radiant life: 9000 h,
 - Disinfectable area: 12 - 15 m²,
 - Lamp rotation (exposure angle can be adjusted): 270°,
 - Anti-schock safety class: I,
 - Case type: IP 20,
 - Operation type: continuous,
 - Chamber's dimensions: 960 x 85 x 135 mm,
 - Chamber's weight: 3 kg,
 - The extension's length: 1800 mm,
- It has a valid Technical Passport,
- Guarantee: 12 months,

