



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NBV 2x30 P germicidal lamp (New)

- A brand new product,
- Made in Poland,
- Direct-acting germicidal lamp NBV 2x30 P: without timer,
- Irreversibly deactivates viruses and microorganisms,
- Equipped with a reflector made of reflective aluminum, which guarantees the highest quality of the device,
- **The use of the device:**
 - Operating rooms,
 - Treatment and bandaging rooms,
 - Patient rooms,
 - Clinics,
 - Laboratories,
 - Pharmacies,
 - Pharmaceutical industry,
 - Food industry,
 - Cosmetics industry,
- **Advantages of NBV 2x30 P germicidal lamp:**
 - Mounting type: M – mobile, on a tripod,
 - The purpose of the product is to provide a level of microbiological cleanliness of the interiors,
 - These devices emit UV-C radiation with a wavelength of 253.7 nm,
 - They irreversibly deactivate bacteria, viruses, molds, fungi and all other microorganisms,
 - The germicidal lamp effectively disinfects the air and surface,
 - The lamp reflector is made of high quality reflective aluminum,
 - The products are designed for the use in operating rooms, treatment rooms, bandaging rooms, patient rooms, clinics, laboratories, pharmacies, pharmaceutical, food and cosmetic industries,
 - The effectiveness of germicidal lamps is supported by tests,
 - The products have all the necessary certificates,
- **Technical data:**
 - The lamp is placed on a mobile tripod,
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 66W,
 - UV-C emitting element: 2 x 30W,
 - UV-C irradiance at 1 m distance: 3.6 W / m²,
 - Radiant life: 9000 h,
 - Disinfectable area: 18 - 22 m²,
 - Lamp rotation (exposure angle can be adjusted): 270°,
 - Anti-schock safety class: I,
 - Case type: IP 20,
 - Operation type: continuous,
 - Chamber's dimensions: 960 x 85 x 145 mm,
 - Chamber's weight: 5 kg,
 - The extension's length: 1800 mm,
- It has a valid Technical Passport,
- Guarantee: 12 months,

