



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **NBV 2x30 PL germicidal lamp (New)**

- A brand new product,
- Made in Poland,
- Direct-acting germicidal lamp NBV 2x30 PL with timer,
- The device has a counter with a microprocessor, without a display, with an acoustic indication when the lamp is replaced,
- It has a CL-02 reader to read the current status of the counter,
- **The use of the device:**
  - Medical institutions,
  - Hospitals,
  - Pharmaceutical industry,
  - Food industry,
  - Cosmetics industry,
- **Advantages of NBV 2x30 PL germicidal lamp:**
  - Counter type: L - digital counter with a microprocessor, without a display. With the use of sound, this model signals the need to replace the radiator,
  - Mounting method: M - mobile on a tripod,
  - They are designed to raise and maintain the level of microbiological purity of rooms,
  - These devices emit UV-C radiation with a wavelength of 253.7 nm,
  - They irreversibly deactivate bacteria, viruses, molds, fungi and all other microorganisms,
  - The germicidal lamp effectively disinfects the air and surface,
  - Reflector made of reflective aluminum with very high quality,
  - Designed for use in operating rooms, surgery rooms, bandaging rooms, patient rooms, clinics, laboratories, pharmacies, pharmaceutical, food and cosmetic industries,
  - Effectiveness of germicidal lamps confirmed by tests,
  - All required certificates,
- **Technical data:**
  - The lamp is placed on a mobile tripod,
  - Supply voltage: 230 V 50 Hz,
  - Power consumption: 66 W,
  - UV-C emitting element: 2x30 W,
  - UV-C irradiance at 1 m: 3.6 W / m<sup>2</sup>,
  - Radiant life: 9000 h,
  - Disinfectable area: 18 - 22 m<sup>2</sup>,
  - Lamp rotation (exposure angle can be adjusted): 270°,
  - Anti-schock safety class: I,
  - Case type: IP 20,
  - Operation type: continuous,
  - Chamber's dimensions: 960 x 85 x 145 mm,
  - Chamber's weight: 5 kg,
  - The extension's length: 1800 mm,
  - Counter type: without a display,
- It has a valid Technical Passport,
- Guarantee: 12 months,

