



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Steriliser MEDICAL MAXI MEDIQUARTZ 138-7 (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in Italy,
- Power supply: 230 V,
- Frequency: 50 - 60 Hz,
- Power: 170 W,
- Temperature: 230°C,
- Container: diameter: 6,5 cm depth: 5 cm,
- Thermometer,
- Weight: 2,5 kg,
- Dimensions: 13,5 x 15,5 x 13,5 cm,
- Technical Passport,
- Warranty: 30-days

