



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Videoprocessor and Light Source OLYMPUS OTV-SI (Reconditioned)**

Used, perfect technical condition,

Japanese production,

Power supply: 230 V,

Frequency: 50-60 Hz,

Power: 250 VA,

Adapted for flexible and rigid endoscopes,

Foiled, tactile control panel,

Lighting: halogen bulb 15V 150W 2pcs.

Possibility to switch the bulb (in case of burning) from A to B,

Manual brightness adjustment,

Included:

OLYMPUS keyboard,

BNC cable,

S-VIDEO cable,

Power supply cable,,

Size: 29,5 x 44 x 17 cm,

The manual is in English (pdf format),

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan.

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

