



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



AED PROGETTI RESCUE SAM (New)

- New,
- Made in Italy,
- Portable with handle,
- Biphasic,
- Semiautomatic public use defibrillator guides the rescuer through the cardiopulmonary resuscitation operations. Thanks to its sophisticated analysis algorithm, fully developed in Italy, this important device knows when the electric shock is needed,
- Compact, ergonomic and lightweight, RESCUE SAM is a life-saving device, indispensable to make your environment a cardio-protected site,
- **Battery:** 12 V, DC, 4,2 Ah, Lithium-ion battery single use,
- **Capacity:** 200 shocks or 4 h continuous work,
- **Voice commands,**
- **Charging time:** 8 sec with new battery charging time can extend when temp is below 10°C,
- **Energy:** 200 J delivered in 50 ohm charging,
- **Easy to use:** by single button SHOCK,
- **Everyday auto test,**
- **Comes with:**
 - AED,
 - Adult paddles,
 - Battery,
- **Additional equipment (extra paid):**
 - Pediatric paddles,
 - Carrying case,
 - Software with USB cable,
 - Wall-mount,
 - Wall locker,
- Weight: 2,2 kg,
- User manual,
- Dimensions: 30 x 28 x 9 cm,
- Technical Passport,
- 24-month warranty

