



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Patient transfer trolley LUNA 02 (New)

New,
Made in Poland,
LUNA-02 trolley is intended to transport a patient in supine position within the hospital or other medical facilities,
The trolley is made of powder-coated steel profiles,
The height of the trolley is adjusted hydraulically,
The top of the trolley consists of HPL plate transparent to X-rays and polyurethane foam mattress covered with MEDITAP fabric,
In standard the trolley is equipped with:
Folding side rails,
Adjustable headrest with a gas spring,
Bumper disks,
As additional options (at extra charge) available:
Central caster blocking,
IV hanger,
Oxygen bottle holder,
Basket for personal belongings,
Patient restraint straps,
Mattress with handles for patient transfer,
Guide for X-ray cassette,
Frame covered with antibacterial paint in white,
Technical data:
Mattress length/total: 195/218 cm,
Mattress width /total: 65 /78 cm,
Height: 47-86 cm,
Headrest angle adjustment: 0° to 80°,
Maximum load: 200 kg,
Weight: 80 kg,
It has a valid Technical Passport,

Warranty: 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

