



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Patient transfer trolley LUNA 03 (New)

- New
- Made in Poland
- TC-11 intended for patient transfer in horizontal position
- Frame made from powder coated steel
- Hydraulic height adjustment
- Mattress made from polyurethane covered with MEDITAP material
- Smooth Trendelenburg due to gas pump
- Standard Equipment:
 - side rails
 - adjustable headrest (gas spring)
 - bumpers
- **Additional Equipment (paid):**
 - wheels with central break
 - IV Pole
 - handle for gas cylinder
 - safety belts
 - basket for accessories
 - X-ray
 - mattress with handles
 - frame covered with white antibacterial coat
- Technical Data:
 - lenght mattress/total : 195/218 cm,
 - width mattress/total: 65 /78 cm,
 - Height: 47-86 cm,
 - headrest angle: 0° do 80°,
 - Trendelenburg: 0° do 12°,
 - Max lifting capacity: 200 kg,
- Technical Passport
- 12 months warranty

