



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Patient transport trolley LUNA 04 (New)**

New,  
Of Polish manufacture,  
LUNA 4 trolley is intended to transport a patient in supine position within the hospital or other medical facilities,  
The trolley is made of powder-coated steel profiles,  
The height of the trolley is adjusted hydraulically,  
The top of the trolley consists of an HPL board transparent to X-rays and a polyurethane foam mattress covered with MEDITAP material,  
The patient transport trolley has smooth regulation of Trendelenburg tilt with a possibility of blocking with a gas spring,  
As a standard the trolley is equipped with:  
Folding side rails,  
Adjustable headrest with a gas spring,  
Bumper disks,  
As additional options (at extra charge) are available:  
Central caster blocking,  
IV hanger,  
Oxygen bottle holder,  
Basket for personal belongings,  
Patient restraint straps,  
Mattress with handles for patient transfer,  
Guide for X-ray cassette,  
Frame covered with antibacterial paint in white,  
Technical data:  
Mattress length/total: 195/218 cm,  
Mattress width /total: 65 /78 cm,  
Height: 47-86 cm,  
Headrest angle adjustment: 0° to 80°,  
Trendelenburg position: 0° to 12°,  
Maximum load: 200 kg,  
Kerb weight: 83 kg,  
Holds a valid Technical Passport,  
Warranty: 12 months,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

