



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Patient transfer trolley LUNA 05 (New)

New,
Of Polish manufacture,
LUNA 5 trolley is intended to transport a patient in supine position within the hospital or other medical facilities,
The trolley is made of powder-coated steel profiles,
The height of the trolley is adjusted hydraulically,
The top of the trolley consists of HPL board transparent for X-rays and polyurethane foam mattress covered with MEDITAP fabric
The patient transport trolley has smooth regulation of Trendelenburg and Anti-Trendelenburg tilt with a possibility of blocking with a gas spring,
In standard the trolley is equipped with:
Folding side rails,
Adjustable headrest with a gas spring,
X-ray cassette guides,
Bumper disks,
Shelf for patient's personal belongings,
As additional options (at extra charge) available:
Central caster lock,
IV hanger,
Oxygen bottle holder,
Patient restraint straps,
Mattress with handles for patient transfer,
Frame covered with antibacterial paint in white,
Technical data:
Mattress length/total: 195/218 cm,
Mattress width /total: 65 /78 cm,
Height: 57-93 cm,
Headrest angle adjustment: 0° to 80°,
Trendelenburg position: 0° to 15°,
Anti-Trendelenburg position: -12° to 0°,
Maximum load: 200 kg,
Kerb weight: 85 kg,
Holds a valid Technical Passport,
12 month warranty,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

