



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Probe for Sonotherapy GS-1/1 (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in Poland,
- Operating frequency:
 - 1 MHz or 3,5 MHz,
- Max power: 2,5 W,
- Max power density: 2,5 W/cm²,
- Range: 1 cm²,
- Dimensions: 18 x 3,5 x 6,5 cm,
- **Compatible with:**
 - ASTAR SONARIS, SONARIS S,
- **Comes with:**
 - probe handle,
- Technical Passport,
- Warranty: 30-days

