



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Rhinolaryngoscope OLYMPUS ENF P3 (Reconditioned)**

- Condition: Reconditioned - Very Good,
- Made in Japan,
- 4 black spots,
- Direction of View: Forward,
- Outer diameter insertion tube: 3,6 mm,
- Working length: 260 mm,
- Max Angulation Up/Down: 130°/130°
- Field of View: 85°,
- Original suitcase: 42 x 32 x 9 cm,
- Technical Passport,
- Warranty: 30-days

