



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Needle Heyinovo WS 2418 PN 2305 (disposable) (Ex-Display)

- Condition: New: Ex-Display,
- Diameter: 2,4 mm,
- Length: 180 cm,
- Needle diameter: 23G,
- Needle length: 5 mm.

