



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Mammography Chair AKRUS AK 5010 MBS (New)

- New,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Application:
 - Mammography,
 - Biopsy,
 - Tomosynthesis,
 - Stereotactic,
- Adjustable back support,
- Height adjustable head rest,
- Drop down shoulder segments and arm rests,
- Safety switch to lock height adjustment during biopsy procedure,
- Foot controls for castors on both sides,
- Fold out foot support,
- Infinitely adjustable,
- Biopsy position,
- Shock position,
- Safe and comfortable support in "park bench" position,
- Electrical height adjustment: 680 to 980 (mm),
- Dimensions:
 - Length wheel base (backrest up) 760 mm,
 - Width wheel base 580 mm,
 - Width chair incl. Armrest 760 mm,
 - Width seat cushion 620 mm,
 - Height chair back rest upright 1440 mm,
 - Max. length chair, back rest horizontal 1680 mm,
 - Wheel diameter 125 mm,
 - Max. patient weight 135/250 kg,
 - Chair weight (depending on options) 95 kg,
- Additional Equipment:
 - Back stabilizer,
 - Shifting head rest,
 - Elastic cushions,
 - Accumulator,
- Technical Data,
- 12-month warranty

