



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Doctor Stool AKRUS AK 445 - Pneumatic height adjustment by hand release (New)**

- New
- Made in Germany
- Round seat
- Backrest
- Multi articulated arm rest, short. or extended
- Pneumatic height adjustment by hand release
- On castors,
- Dimensions:
  - Base diameter: 600 mm,
  - Castors: 50 mm,
  - Adjustable height: 530 - 720 mm,
  - Seat diameter: 420 mm,
  - Seat thickness: 80 mm,
  - Max load: 135 kg,
- Optional: foot rest, chrome-oval
- Technical Passport
- 12 months warranty

