



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Doctor Stool AKRUS AK 445 - Pneumatic height adjustment by hand release (New)

- New
- Made in Germany
- Round seat
- Backrest
- Multi articulated arm rest, short. or extended
- Pneumatic height adjustment by hand release
- On castors,
- Dimensions:
 - Base diameter: 600 mm,
 - Castors: 50 mm,
 - Adjustable height: 530 - 720 mm,
 - Seat diameter: 420 mm,
 - Seat thickness: 80 mm,
 - Max load: 135 kg,
- Optional: foot rest, chrome-oval
- Technical Passport
- 12 months warranty

