



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Operator's seat with backrest AKRUS AK 445-football height adjustment (New)

New,
German production,
AK 445 is a comfortable operator's seat with a round seat,
Built-in backrest,
Manual height adjustment,
Mobile, five-fingered base with wheels,
Overall dimensions:
Diameter of the base 600 mm,
Wheels with a diameter of 50 mm,
Height adjustment 530 - 720 mm,
Diameter of seat 420 mm,
The thickness of the seat 80 mm,
Max. the weight of the person sitting 135 kg,
Has a current Technical Passport issued,
Warranty 12 months

