



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Surgeon Chair AKRUS AK 480 (New)**

- New
- Made in Germany
- Power supply: 230V or battery
- Frequency: 50/60Hz,
- Surgeon chair, electric height adjustable, offering a number of sensible and ergonomically valuable options.
- All standard functions or additional options have been designed to support a fatigue free seating posture of the surgeon preventing any long term strains due to improper seating positions.
- Seat cushion tilt by approx. 12°
- The height adjustable back rest is spring loaded, but also lockable in position.
- Heavy duty dual castors , lockable by one foot heel kick
- Electrically height adjustable by heel kicked foot switches
- Optional Equipment:
  - The upper part of the chair moves 100 mm fwd – bwd, conveniently controlled by a foot release lever
  - The standard arm rest features height adjustment and rotation.
  - The height adjustable arm rests pivot and tilt to any desired position on a ball joint
  - The height adjustable arm rests pivot and tilt to any desired position on a ball joint with extended arm
  - The angle of the moving lateral wings is adjustable.
  - Additional battery
- Technical Data:
  - Base diameter: 620 x 680 mm,
  - Adjustable height: 550 x 720 mm,
  - Seat: 420 x 420 mm,
  - Seat thickness: 80 mm,
  - Adjustable armrest: 140-290 mm,
  - Castor diameter: 75 mm,
  - Weight: 55 kg,
  - Max load: 250 kg,
- Technical Passport
- 12 months warranty

