



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Instrument Delivery Table AKRUS AK 106 vario (New)

- New,
- Made in Germany,
- Power supply 230 V,
- Frequency 50/60 Hz,
- Affordable Instrument delivery table,
- Electric infinitely variable height adjustment,
- Assembly kit designed to position lift at center or off center position,
- Durable, scratch resistant laminated table top,
- 4 stands or 4 castors,
- 2 IEC power outlets,
- Technical Data:
 - Base: 400 x 600 mm,
 - Adjustable height: 590 – 890 mm,
 - Max load: 50 kg,
 - Weight 39 kg,
- Countertops:
 - Small top 50 x 60 mm,
 - Large top 60 x 100 mm,
 - V shape 60 x 100 mm,
- Technical Passport,
- 24-month warranty

