



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Examination Chair AKRUS SC 5010 HS (New)

New,
German production,
Power supply 230V,
Frequency 50/60Hz,
Mobile operating table designed for e.g:
ENT surgery,
Maxillofacial surgery
Plastic surgery,
Split seat - movable leg support
Electrically adjustable headrest, backrest, seat and leg rest,
Each element individually adjustable,
Adjustable height - possibility of extremely low setting 64-100 cm,
Seat level adjustment,
Possibility to mount various ergonomic headrests,
Mobile - on 125mm castors with brakes on both sides,
Possibility of battery operation,
Large selection of additional accessories,
Dimensions:
Length of wheel base 760 mm,
Width of the base of the wheels 580 mm,
Table width plus arm 700 mm,
Seat width 590 mm,
Height of backrest with upright backrest 1700 mm,
Max. Length, backrest horizontal 1750 mm,
Max. Patient weight 200/300 kg,
Weight (depending on options) 114 kg,
Optional equipment (at extra charge):
Standard headrest,
Multi articulating headrest,
Flat headrest,
Articulated arm for additional headrest +30/-50°,
Screen hanger,
Ventilated screen holder,
Armrest, adjustable in height, pivoting, on ball head,
Arm rest, not adjustable,
IV-Pole,
Side rail,
Foot switch,
Foot joystick 5 segments,
Battery,
Holds the latest Technical Passport,
Warranty 24 months,
Finance options (Poland only): Instalments, Leasing, Loan

