



PRAISTON Sp. z o.o.
ul. Górska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Autoclave HYDRA 100 (New)

- New,
- Made in Italy,
- Class N,
- Features:
 - The autoclave HYDRA 100 is equipped by a heating element placed under the sterilization chamber which allows the heating to 125°C of the distilled water. The pressure that is created inside the sterilization chamber, in combination with the temperature, allow the complete asepsis on the instruments,
 - The sterilizer is equipped with several sensors that monitor the reaching of the temperature, which is preset and can not be changed to prevent tampering,
 - The machine is equipped with several safety systems, in particular with a thermostat that, in the event of an excessive temperature, stops the sterilization cycle. It is also equipped with a safety valve with TÜV certification (PED IV), the same installed on the class B autoclaves,
 - With little maintenance and the replacement of a few parts subject to wear this autoclave can continue to sterilize for more than 15 years,
 - Available in version 9 liters and is compliant with Directive 93/42/EEC on small steam sterilizers,
 - With its reliability, speed, and less space is a good alternative to class B autoclaves for the sterilization of unwrapped solid instruments,
- Technical Data:
 - Class N,
 - Power supply: 220/230 V,
 - Frequency: 50 - 60 Hz,
 - Power: 1500 W,
 - Temperature: 125°C,
 - Operating pressure: 1,3 - 1,6 bar,
 - Dimensions: L 250, H 310, W 385 mm,
 - Chamber dimensions: ø220 - H 235 mm,
 - Distilled water,
 - Chamber capacity: 9 litres,
 - Time per cycle: 15 minutes,
 - Weight:
 - Net: 13,3 kg,
 - Gross: 17,5 kg,
- Technical Passport,
- 12-month warranty

