



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical equipment table TALIE 4 (New)

New,

Made in Poland,

It has a stable steel construction covered with powder paint,

Standard features include an upper metal shelf with rims and holders for electrodes and suction cups, as well as casters, 2 of which have a brake,

As additional options are available (at an extra charge):

Additional metal shelf,

Shelf with drawer and stainless steel bowl,

Gel holder,

Paper hanger,

Technical parameters:

Width: 53 cm,

Height: 89 cm,

Depth: 42 cm,

Shelf dimensions (W/D): 36x 30 cm,

Dead weight: 14 kg,

Maximum load per shelf: 5 kg,

Warranty: 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

