



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Respiratory gas monitor DRAGER PM 8050 (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50-60 Hz,
- Power: 80 VA,
- LCD Display with brightness adjustment,
- Foiled control panel,
- Menu in English,
- Readable and easy to use menu,
- **Comes with following modules:**
 - SPO2 - not included (equipping extra paid),
 - TEMP - not included (equipping extra paid),
 - O2 - not included (equipping extra paid),
 - SPIRO - not included (equipping extra paid),
- Gas measurements: DES, ENF, HAL, ISO, SEV, as well as: O₂ i N₂O, CO₂,
- DRAGERWaterLock i new oxygen sensor,
- Inputs: RS 232 C, Monitorbus, analog output
- Acoustic and visual alerts,
- Monitor dimensions: 42,5 x 33 x 19,5 cm,
- User manual in English (pdf),
- Technical Passport,
- Warranty: 30-days

