



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Headrest FAMED WS 46-5 (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in Poland,
- Used for supporting a patient's head in a supine position and lateral position,
- The headrest is height, angle and longitud adjustable,
- Max load: 15 kg,
- Headrest size: 50 x 27,5 x 5 cm,
- Technical passport,
- Warranty: 30-days

