



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Arthroscopy set AESCULAP (Reconditioned)

Used, very good condition,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
The set includes:
NEOVO Model X-19AV monitor,
SONY DVD Writer Model RDR-AT200,
DVD recorder with hard drive capacity: 250 GB,
Record format: -R (Video) , -R (VR) , +R , -RW , +RW , -R (Video) DL (two layer) , -R (VR) DL (two layer) , +R DL (two layer) , -RW DL (two layer) , +RW DL (two layer),
Analogue terrestrial TV tuner: PAL-B/G, D/K, I, SECAM-L, SECAM-D/K,
Entrances:
Composite video signal 1 x front,
S-Video 1 x front,
Audio input (L/P) 1 x front,
i.LINK™ (DV) input 1 x front,
RGB recording input (LINE1 SCART),
USB type A 1 x front (HDD camera, digital camera, flash memory, hard disk),
USB type B compatible with PictBridge 1 x front,
Aesculap ProMIS Line Fluid 2 pump Model: PG-030,
Max. Flushing and suction power 500mmHg,
Max. flow 1000 ml/min,
Aesculap Flow 25 CO2 Inflator Model: PG-060,
Aesculap CCD-CAMERA Camera Control Unit Model PV-280,
CCD module,
Endoscopic camera processor adapted to rigid endoscopes,
White balance control,
The "small image range" function,
Foiled control panel,
Aesculap camera Model: PV121S,
Light source AESCULAP Model VIDEOLIGHT 400 OP 926:
Used, very good condition,
German production,
Adapted for rigid endoscopes,
Halogen light,
Power supply: 230 V,
Frequency 50-60 Hz,
Power: 600 VA
Light bulbs:
250 W metal halide lamps - proper lighting,
Halogen 250 W - emergency lighting,
Colour temperature: 54000K,
Six-stage intensity control,
Equipped with a carrying handle,
Foiled control panel,
Size: 35 x 30 x 32.5 cm,
Powershaver SL Karl Storz, (without accessories, aftermarket purchase possible),
SPOINA endoscopic trolley,
Used, very good condition,
Polish production,
Mobile on wheels 2 with brake,
It has 4 shelves,
A camera head holder,
Lockable locker,
Space for gas cylinders,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

