



Rehabilitation hospital - bed MALVESTIO 3LE100 (Ex-Display)

New, exhibition,
Italian production (2009),
Bed for adult patients staying in hospitals during diagnostics, monitoring and treatment,
Four-section hospital bed with a fixed height,
Laid segments (galvanized grille),
Technical data:
Dimensions 2160 x 930 mm,
The angle of elevation of the 75 ° backrest segment is adjusted by means of a gas actuator,
Thigh segment angle 22 ° adjustable by means of a gas cylinder,
Wheel diameter 125 mm,
Tripod for 4-arm drip,
GOOD anti-bedsore mattress from JANPOL made of polyurethane foam with a density of 25 kg / m³, height 120 mm with a pair - permeable - watertight,
Maximum load 170 kg,
Mobile on wheels with brakes
Has a current Technical Passport issued,
1-month warranty

