



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DIRUI H-500 Urine analyzer (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in PRC,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 40 VA,
- DIRUI H-500 is a semi-automatic strip reader for urinalysis,
- Features:
 - Large 5.7" LCD touch screen, convenient for users to operate,
 - Software in polish,
 - Built-in thermal printer (57mm wide paper),
 - Environment friendly system for collecting and disposing used strips,
 - Barcode reader,
- Specifications:
 - Cold light source with 4 wavelengths ensures accuracy of results,
 - Measuring system: diode array, monochromatic reflected light measurement,
 - 514 samples per hour,
 - 10 or 11 test parameters,
 - Strips used: H10 and H11,
 - Memory of 2000 results, convenient to review,
- Included:
 - Power cord,
- Overall dimensions: 37 x 40 x 30 cm,
- Weight: 7 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

