



Blood cell mixer Ortho BIOVUE SYSTEM (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50-60 Hz,
- Power: 25 W,
- Foiled control panel with digital display,
- Centrifuging: full cycle - 5 mins,
- Electronic cover lock,
- Inductive brushless engine designed for a long-term use without maintenance,
- Chamber diameter: 24,5 cm,
- Rotor with 12 cassette sockets,
- Dimensions: 32,5 x 37,5 x 24 cm,
- Technical passport,
- Warranty: 30-days

