



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RODENSTOCK PRO 1000 Evo Ophthalmic unit (New)

- New,
- Made in Germany,
- The PRO 1000 EVO ophthalmic unit has been equipped with new features:
 - Electric height adjustment of the unit's table,
 - Electro-magnetic table brake,
- Features and advantages of the RODENSTOCK PRO 1000 Evo Unit:
 - Robust design,
 - Comfortable, adjustable seat,
 - Timeless, modern design,
 - Column with top LED lighting,
 - Internal RGB LED lighting,
 - Movable table for 2 devices,
 - Mechanical locking of the table,
 - Adjustable phoropter arm,
 - Projector stand,
 - Possible additional cabinet,
- Technical specifications:
 - Rotating and sliding table, for two devices,
 - Locking of the table: magnetic,
 - Table adjustment: +/- 10 cm
 - Max load capacity (table): 45kg (25/20kg),
 - Patient seat: ER 500,
 - Patient seat elevation: ± 20 cm,
 - Patient seat tilt: up to 33° ,
 - Total weight (without chair): approx. 195/215 kg,
 - Max load capacity (chair): 120 kg,
 - Projector column height: 176 cm/ ± 1.5 cm,
 - The unit is equipped with the following outlets:
 - 2 power sockets (230 V),
 - 1x RS-232, 1x 6-pin (6 V/12 V),
 - 1x LAN,
 - 1x USB 2.0,
 - Reading light: 12 V white LED,
 - Interior light: RGB LED 12 V,
 - Power consumption: 280 VA,
 - Power supply: 110/115 V AC/230 V AC $\pm 10\%$,
- Has a current Technical Passport issued,
- Warranty 12 months,
- Valid technical passport issued,

