



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ECG FUKUDA CARDISUNY ALPHA 1000 (Reconditioned)

Used, good technical condition,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 30 VA,
portable
3 channel,
Built-in printer paper width 60mm,
Write speed: 25/50 mm / s,
Included:
FUKUDA KP-313 ECG cable 10-lead,
Overall dimensions: 27.5 x 20.5 x 7 cm,
Has a current Technical Passport issued,
1-month warranty

