



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Ultrasonic probe BARD DYMAX (Reconditioned)

Used, perfect technical condition,
Production U.S.A (2006),
Frequency of work:
9 MHz 53 mmFL,
Dimensions: 13.5 x 4 cm,
Meant for:
Vascular USG BARD SITE RITE IV,
Has a current Technical Passport issued,
1-month warranty

