



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Laser LASERSCOPE STONELIGHT (Reconditioned)**

Used, perfect technical condition,

U.S.A production,

Date of production: 2009,

Power supply: 230 V,

Frequency: 50-60 Hz,

Power: 10 A,

StoneLight laser - is a versatile device, applicable in various applications, in particular in urology. It is intended for use in procedures such as cutting, ablation and vaporization. Depending on the specifics of the procedure, you can use contactless methods, operate with the help of an endoscope on the open field. It is a breakthrough solution in the treatment of kidney stones, in the removal of malignant lesions or deposits.,

Laser type: Holmowo-Yag rod (THC-YAG),

Light length: 2100 nm,

Impulse duration: 350 ms,

Frequency of pulses: 5,8,10,12 PPS (impulses per second),

Power range: 1-15 W,

Pulse energy: adjustable: max 1.5 J,

Exposition: Single pulse or continuous operation,

Target beam: laser diode, red, 635 nm,

Built-in display with foil panel,

Retractable fiber optic support,

The convenient sequence of connections allows you to choose the optimal fiber length to work: 150, 200, 320, 550  $\mu$ m. The ability to sterilize fibers.

Included:

The main LASERSCOPE control unit,

3 pieces of protective glasses with a case,

Footbridge,

Remote lock attachment,

2 keys for switching on the laser,

Power cord,

Mobile on wheels with brakes,

Weight: 59 kg,

Dimensions: 36 x 53 x 93 cm,

User manual in Polish, English,

Has a current Technical Passport issued,

1-month warranty

