



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RODENSTOCK RSL 4300 Slit lamp (New)

- New,
- Made in Germany,
- Advantages and features of the RSL 4300 Slit Lamp:
 - Zeiss type lamp (lower),
 - Halogen or LED illumination,
 - Stereomicroscope with three magnification ranges,
 - Converging binocular microscope,
 - Exceptional optical quality,
 - Special FLUO cartridge for optimal observation of fluorescent effect,
 - 4 filters in the projector: heat-absorbing, green, cobalt green, UV,
 - Slit length: 0 -12 mm,
 - Built-in yellow filter,
 - Complete with chin, transformer + accessories,
 - BLACK series,
- Specifications:
 - Microscope:
 - Version with three magnifications,
 - 12.5x magnification eyepiece,
 - Total magnification factor: x10, x16, x25 ,
 - Pupil distance: 49 to 77 mm,
 - Focal distance: 107 mm,
 - Field of view: 22, 14, 8.5mm,
 - Convergence angle: 13°,
 - Slit illumination:
 - Slit width: 0 - 14mm stepless,
 - Slot length: 1 - 14mm,
 - Hole diameters: 0.2, 1 mm square, 2, 3, 5, 9, 14 mm,
 - Filters: heat absorption filter, blue filter, red-free filter,
 - Horizontal focus adjustment: 12 mm,
 - Light source: LED / 12V,
 - Additional features:
 - Additional diffuser,
 - Yellow filter integrated into the microscope,
 - TAT-100R applanation tonometer,
- Dimension and weight:
 - Dimension (WDH): 350 x 432 x 604 mm,
 - Weight: 19 kg,
- Power requirements:
 - Power supply: (100-240 V input) +/- 10% multi-plug complying with EN 60601-1, EN 61000-6-2, EN 61000-6-3,
 - Power supply: 12V DC / 2.5 A IEC / EN 60 6001,
- It has a valid Technical Passport issued,
- Warranty 12 months,

