



RODENSTOCK PERSMART Perimeter (New)

- New,
- Made in Germany,
- Product advantages and features:
 - Dedicated to the examination of drivers,
 - Innovative solutions,
 - Static measurement,
 - LED technology,
 - Simple and fast operation,
 - Small and compact housing,
 - Electric chin,
 - Modern and precise technology,
 - Digital eye tracking,
 - White and white stimulus,
- Specifications:
 - Type of stimulus bowl: Hemispherical part,
 - Bowl radius: 300 mm,
 - Measuring range: 50° horizontally (+30° with movable attachment point), 30° vertically (+10° with movable attachment point),
 - Chin adjustment: Electrically controlled,
 - User interface: Multilingual,
 - STANDARD TEST FIELDS: Full, Glaucoma, Central 30°, Central 22°, Spot 10°.
 - SPECIALIZED TESTS: Bi-driving, industrial medicine, augmented industrial medicine,
 - SOURCE OF BODDER: LED rear projection,
 - SOURCE COLOR: White,
 - SOURCE SIZE: Goldmann size III (0.43°),
 - BACKGROUND DENSITY: 0.01 cd/m² (0.03 asb) to 333.33 cd/m² (1000 asb) in 3 dB or 1 dB increments
 - BACKGROUND ILLUMINATION: Automatic level control 10 cd / m² (32asb),
 - EXPOSURE TIME: Adjustable from 0.1 to 9.9 sec,
 - FIXATION CONTROL: Fixation control: Heijl-Krakau blind spot monitoring, digital eye tracking,
 - DIMENSIONS AND ELECTRICAL REQUIREMENTS:
 - Dimensions (WxDxH): 548 x 382 x 450 mm
 - Weight: Approximately 9 kg,
 - Power supply: 115 to 230 V AC,
 - Frequency 50/60 Hz,
 - Power consumption: Max 35 VA,
- Has a current Technical Passport issued,
- Warranty 12 months,

