



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Ankle brace BAUERFEIND AirLoc (Reconditioned)**

Used, very good technical condition,  
German production,  
indications

Acute trauma to the joints and ligaments of the joint within the ankle,

Chronic instability,

Postoperative rehabilitation,

Traumatic prevention,

Features

AirLoc is a stabilizing brace with an innovative flexible shell-shaped coating. This ensures a perfect fit to the shape of the leg and the size of the swelling,

Plastic, anatomically shaped cups ensure excellent adhesion to the hocks,

The quadruple velcro fastening system allows individual adjustment which guarantees the highest comfort of use,

Air-filled panels allow easy adjustment of the orthosis to the size of the swelling,

The AirLoc system ensures a high degree of stabilization of the ankle,

The slender cups work well when wearing an orthotics in shoes,

Action

It supports sensory-motor functions while preventing the ankle from twisting,

Stabilizes the lateral capsule ligaments of the upper ankle, especially if it is worn inside the shoe,

Has a current Technical Passport issued,

Warranty for 1 month

