



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Ankle stabilizer with flexible rails NEOX S-03 (Ex-Display)

Used, very good technical condition,
Korean production,
Power supply: 230 V, through an 18V power supply,
Frequency: 50-60 Hz,
Power: 1.5 A,
7 '(inch) high resolution color LCD display with adjustable angle of inclination, and 360 ° rotation,
Memory: 450 hours of recording,
Possibility to connect to a computer,
Built-in connector: RS-232,
Built-in thermal printer paper width: 90 x 150 mm,
FHR - range 30-240 bpm,
Write speed: 10, 20, 30 mm / minute,
Auto Print: OFF / 10/20/30/40/50/60 minutes,
3 channel screen,
The set includes:
1 x TOCO head,
1 x US head,
Systolic button of the patient,
Power cord + power supply,
5 packages of paper,
Camera dimensions: 320 x 280 x 105/230 mm,
Has a current Technical Passport issued,
1-month warranty

