



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Universal scales (baby-adults) TANITA 1582 (Reconditioned)

Used, very good technical condition,

Japanese production,

Power supply: 6 1.5 V AA batteries,

Digital display with rotation display result display,

Load max. 136 kg,

Accuracy:

30lb / 14 kg 1oz / 20g,

300lb / 136kg, 0-250lb / 0,2lb 250-300lb 0,5lb, 0-100kg / 100g, 100-136kg / 200g,

Unit of measure: lb / kg,

It has two functions: the neonatal weight or after dismantling the pan serves as a floor scales weighing up to 136 kg,

Dimensions of the base: 320 x 340 x 50 mm

The dimensions of the pan: 565 x 275 mm,

Weight: 3.7 kg,

Has a current Technical Passport issued,

1-month warranty

