



### **Tabletop peristaltic pump CARE FUSION ALARIS SE (Reconditioned)**

Used technical condition is perfect  
made in Poland (2011)

Power supply: 230 V and through the built-in NiCd battery (1800mAh),

Frequency: 50-60 Hz

Power: 75 VA,

each of the Alaris SE pumps are equipped with a regulator flow AccuSlide ® . AccuSlide minimizes the risk of free flow due to the built-in ogranicznikowi flow which closes automatically when the headset,

Specification:

flow Range: 0.1 to 999.9 ml / HR in 0.1 ml / h (original)

1.1 to: 270.0 cm ml / h in 0.1 ml / HR (secondary),

KVO flow Range: 0.1 to 20.0 ml / h in increments of 0.1 ml / h

volume input: 0.1 to 999.9 ml in 0.1 ml (original) from 1 to 999.9 ml, in a dose of 0.1 ml (secondary)

bulk Density: from 0.0 to 9999.9 ml in 0.1 ml of

Alarms: air in line, battery depleted, channel malfunction, a communication error the computer, disconnected the flow sensor, the shelf life is exceeded, failure of the locking device of the key, opened the latch, there was no significant flow in the stream upuszczanym, down occlusion, penetration of the occlusion is detected when the primary flow,

complete holder for mounting on a tripod, on a bed or other bus, and to move,

RS-232 connector,

Weight: 3.0 kg,

user manual in English (pdf).,

Dimensions: 130 x 220 mm x194,

Has exposed the current Technical Passport.

1 month warranty

