



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### Fluids heater SMITHS HOTLINE LEVEL 1 HL-90 (B) (Reconditioned)

Used, very good technical condition,

U.S.A production,

Power supply: 230 V,

Frequency: 50-60 Hz,

Power: 1.5 A,

Hotline is the only device that heats infused flow through the entire path - from the bag to the patient's injection. Thanks to this, it is the most effective device of this type,

Due to the heating of liquids using the water jacket, there are no so-called 'Hot spots' (as found in dry heaters). The appliance can be used for the safe heating of transfused blood,

Electronic temperature display,

LED work and alarm indicators,

Circulation capacity: 1,4L,

Features of the heater:

Maintains blood and fluids at a temperature of 37 to 42 ° C, providing warm fluids to the patient at a flow rate of up to 5000 ml / h,

Has a holder to hang on the stand,

Consumables available on the Polish market (not included with the device),

Overall dimensions: 21 x 21.5 x 25 cm,

Operating instructions in English (PDF),

Has a current Technical Passport issued,

1-month warranty

