



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LID Floaled Fluorescent light (New)

- New,
- Made in France,
- This LED light illuminates a large area with bright, soft light. The tilting head allows you to easily direct the light and avoid glare,
- The diffused intense illumination of this lamp helps visually impaired people to work more comfortably,
- Specifications:
 - Light intensity: 840 Lumens,
 - Light color: 4500°K* (neutral white),
 - Color rendering index: IRC 95/100,
 - Electricity consumption: 8,4 W,
 - Class 1 medical device,
 - LED lifetime: minimum 50,000 hours,
 - Electrical protection: insulation class II,
 - Standards: IEC 60 601 -1/ -2-41 IEC 60 598-2-25,
 - Power cord length: 3 meters,
 - Available colors : white,
- Dimensions:
 - Lamp head dimension: 10 x 230 mm,
 - Flexible arm: 650 mm, 1000 mm,
 - Base dimension: 69 x 69 x 115 mm,
 - Weight: 1.5 / 2.3 kg,
 - Inter-axis spacing (4 holes in the lamp base): 54.5 mm,
- The standard length of the lamp is 650 mm, for an extra charge you can purchase a lamp with a length of 1000 mm,
- Diameter of the light field :
 - 30 cm distance: 33,000 lux field illumination 15 x 10 cm
 - Distance 50 cm: 13000 lux field illumination 25 x 20 cm,
 - Distance 100 cm: 3500 lux field illumination 40 x 40 cm,
- Additional mounting accessories are required for the lamp , of your choice (at an extra charge):
 - Telescopic stand base 60-100 cm white metal, heavy,
 - Telescopic stand base 60-100 cm white metal, heavy -version with brakes in 2 wheels,
 - Telescopic stand base 90-130 cm metal white, heavy,
 - Telescopic stand base 90-130 cm white metal, heavy - version with brakes in 2 wheels,
 - Ordinary clamp - black,
 - Clamp for medical rails with vertical lamp outlet,
 - Universal clamp,
 - Wall bracket,
 - ARM wall bracket white (2x40cm),
 - ARM wall mount white (1x40cm),
 - Tabletop mount ARM white (1x60cm),
 - Metal base white / gray,
 - Plate with metal pipe in the middle ø16,
 - Plate with metal tube on the edge ø16,
 - Plate on ARM holder,
 - Female socket with adaptor plate Ø 16 mm x 27,
- It has a current Technical Passport issued,
- Warranty 3 years,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

